

Optimizing the intraocular lens formula constant according to intraocular lens diameter

Abdul R El-Khayat¹, Paul Tesha²

¹Department of Ophthalmology, Leicester Royal Infirmary, Leicester, LE1 5WW, UK

²Department of Ophthalmology, Lincoln County Hospital, Lincoln, LN2 5QY, UK

Correspondence to: Abdul R El-Khayat. Department of Ophthalmology, Leicester Royal Infirmary, Leicester, LE1 5WW, UK. a_elkhayat@hotmail.com

Received: 2020-01-15 Accepted: 2021-01-05

Abstract

• **AIM:** To determine whether the different diameters of a specific intraocular lens (IOL) have significantly different optimized SRK/T A constants and whether these new A constants can improve refractive outcomes.

• **METHODS:** Data were collected prospectively from Jan. 2011 to Dec. 2012 on all patients undergoing routine cataract surgery at a district general hospital in the UK. Patients were divided into three groups according to the size of the Akreos AO MI60 IOL used. A constants for the SRK/T formula were optimized according to the size of the IOL. These optimized A constants were then used to select future refractive outcomes.

• **RESULTS:** A total of 2398 cataract operations were performed during the study period of which 1131 met the inclusion criteria. The three optimized A constants for the different sized IOLs were 118.98, 119.13, 119.32. The difference between them was highly significant ($P \leq 0.0001$). Two optimized A constants for three sizes of IOL led to an improvement in refractive outcomes (from 93.4% to 94.6% of refractive outcomes within 1.00 D of predicted spherical equivalent). The optimized A constant for the largest IOL was based on a small number of cases and was not used.

• **CONCLUSION:** Optimizing the A constant for the three distinct sizes of the Bausch & Lomb Akreos MI60 lens lead to three significantly different A constants. In our practice, using two different optimized A constants for three different sized IOLs give the least refractive error, however, using three optimized A constants may give better results with a larger dataset.

• **KEYWORDS:** intraocular lens; formula constant; optimization; cataract surgery

DOI:10.18240/ijo.2021.05.09

Citation: El-Khayat AR, Tesha P. Optimizing the intraocular lens formula constant according to intraocular lens diameter. *Int J Ophthalmol* 2021;14(5):700-703

INTRODUCTION

The position that an intraocular lens (IOL) takes within an eye affects the refractive outcome of that eye^[1-2]. Misalignments of IOLs *in vivo i.e.*, decentration, tilt and axial translation have been shown to impact on the refractive outcome of the eye^[3-4]. Studies have also shown that the configuration of an IOL in the capsular bag depends on three factors: size, shape and rigidity of the IOL^[5]. This suggests that the size of an IOL may affect its configuration/position in the eye and hence the refractive outcome of the eye.

Certain plate haptic IOL designs vary in size in discrete steps across their range of dioptric powers^[6]. This difference in size may affect the IOL position in the eye and hence the accuracy of predictions of postoperative refraction.

Estimating the postoperative IOL position after cataract surgery from preoperative measurements is the main limiting factor for refractive predictability^[7-9]. IOL formulas use constants for specific IOLs to predict the effective lens position postoperatively. These constants are usually refined from refractive outcomes of a large number of cases and hence represent the mean effective position of the specific IOL inside the eye^[10-11]. The optimization of this constant by comparing the predicted post-operative spherical equivalent with the actual spherical equivalent is recommended to improve accuracy of refractive predictions^[12-13].

The main aim of this study was to determine whether optimizing A constants separately for the discrete sizes of IOLs within one design would lead to statistically significant differences in the A constants and whether these can be used to improve refractive outcomes.

SUBJECTS AND METHODS

Ethical Approval Written informed consent for surgery was gained from all participants and the study was conducted according to the principles outlined in the Declaration of Helsinki.

Table 1 Inclusion and exclusion criteria

Inclusion	Exclusion
Cataract operation performed between 1/1/11 - 31/12/12	Complicated surgery
Bausch & Lomb Akreos AO MI60 IOL used	Additional procedures required
	Comorbidities leading to a guarded visual prognosis
	BCVA postop. \leq 6/12
	Postop. astigmatism >3 D
	>2 D difference between actual and predicted spherical equivalents
	IOL not in capsular bag
	SRK/T formula not used

IOL: Intraocular lens; BCVA: Best corrected visual acuity.

Data from the 1st of January 2011 until the 31st of December 2012 from an on-going cataract audit at Lincoln County Hospital were used for this study. These data were collected prospectively on a pro forma and included the following pre-operative and post-operative information: demographics, refraction, visual acuity, keratometry, axial length, A constant used and any surgical complications. Biometry was performed and recorded by different nurses. IOL power calculations were performed according to The Royal College of Ophthalmologists (RCOphth) and The National Institute for Health and Care Excellence (NICE) guidelines (SRK/T formula used for axial lengths of 22 mm and above and Hoffer Q formula used for lengths below 22 mm)^[13-14]. The starting point for the A constant that was used in the SRK/T formula was a constant that had been optimized locally for all sizes of the Akreos AO MI60 lenses collectively. The cataract operations were carried out by the different ophthalmologists who perform this procedure at Lincoln County Hospital. Data relating to choice of IOL implant and complications were recorded by the surgeon. Post-operative refraction was completed by community optometrists at least four weeks after the operation and recorded on a form that was given to the patient upon completion of their surgery. Data collected on all pro-forms were entered into a customised database. Table 1 highlights the inclusion and exclusion criteria for the study.

This study was aimed at highlighting how the size of an IOL affects refractive outcomes, hence only one lens that varies in size was included. The Akreos AO MI60 IOL is the first line IOL used at Lincoln County Hospital and comes in three distinct sizes according to its dioptric power. The different sizes of the IOL are detailed in Table 2.

Analysis of the data consisted firstly of dividing it into three groups according to MI60 IOL size. Then, the difference between predicted post-operative spherical equivalent and actual post-operative spherical equivalent (assessed by community optometrists) was calculated for each IOL size individually. This “error” in prediction was fed back to calculate what, in retrospect, the “ideal” or “optimized” A

Table 2 Different sizes of the Akreos AO MI60 IOL according to dioptric value

Overall diameter (mm)	Dioptric range (D)
11.0	10.0-15.0
10.7	15.5-22.0
10.5	22.5-30.0

constant should have been (*i.e.*, the A constant that would have given the most accurate refractive predictions). This produced three separate A constants for each size of the MI60 IOL. One-way analysis of variance (ANOVA) was used to compare the significance of the difference between the A constants.

Data for errors in predicted post-operative spherical equivalent after A constant optimization were also analysed for the previous 5y. This was done to determine whether using the new A constants would continue to improve refractive predictions. Errors were divided into prediction within 1.5 D, 1.0 D, and 0.5 D of actual post-operative spherical equivalent. Initially, optimization from 2008 until 2010 was carried out for all the MI60 lenses collectively. During 2011 and 2012, two separate A constants were optimized for the medium and small MI60 lens. Optimization was not done for the large MI60 lens as the number of the lenses used was low.

RESULTS

During 2011 and 2012, 2398 cataract operations were performed at Lincoln County Hospital. After applying the inclusion and exclusion criteria, this number was reduced to 1131. Table 3 shows the demographics of the individuals included in the study. One-way ANOVA between mean age, post-operative days to refraction and laterality showed no significant difference between the groups with different sized IOLs.

Table 4 and Figure 1 show the three A constants optimized for IOL size. One-way ANOVA for both mean biometry errors and optimized A constants showed highly significant statistical differences of $P=0.014$ and $P\leq 0.0001$ respectively.

Table 5 and Figure 2 show the accuracy of refractive predictions over the last 5y after using successively optimized A constants. The trend in results showed increasing accuracy

Table 3 Demographics of individuals with respect to size of MI60 IOL

Parameters	Age		Postop. days to refraction		Right eyes	Left eyes
	No.	Mean±SD	No.	Mean±SD		
All	1131	75.57±8.89	783	39.94±25.97	577	554
Large	74	70.67±9.06	57	36.56±14.51	40	34
Medium	590	75.17±8.93	418	40.15±21.04	301	289
Small	467	76.84±8.51	308	40.29±32.80	238	229




Figure 1 SRK/T A constants individually optimized for the 3 different sizes of MI60 lenses.




Figure 2 Accuracy of predicted spherical equivalents with successive optimization of A-constants Results for 2011 and 2012 derived from 2 A-constants optimized according to 2 different sizes of the MI60 lenses.

with successive optimization. This remained the case during 2011 and 2012 when two A constants were optimized and used.

DISCUSSION

This study shows that optimizing the SRK/T A constant for the three distinct sizes of the Bausch & Lomb Akreos MI60 lens leads to three significantly different A constants. Since the

Table 4 Optimized A constants according to IOL size

Diameter of IOL (mm)	Mean biometry error (D)	Optimized A constant
11.0	-0.076	118.98
10.7	-0.023	119.13
10.5	0.075	119.32

optimization process of A constants utilises errors in refractive predictions, in theory, this could mean that using a different sized MI60 lens leads to a different refractive outcome.

This study also shows that using one A constant that is optimized for all MI60 lenses and two constants that are optimized for two sizes of the MI60 lens leads to better refractive outcomes than no optimization. Demonstrating that even better refractive outcomes can be achieved through three different A constants would be ideal, however using three A constants in our practice at Lincoln County Hospital did not lead to improved refractive outcomes and hence is not being used at present.

Lens formula constants published by manufacturers are typically intended for use with contact ultrasound biometry^[15]. Using optical biometry with these constants will lead to outcomes that are more hypermetropic^[12]. Optimized constants for different lenses and methods of biometry are available online and probably provide a better starting point to manufacturer constants^[16]. However, optimization according to local clinical outcomes is still recommended^[15]. Most clinicians who optimize lens constants do so for all eyes collectively^[17]. Optimizing different constants for subsets of patients is more controversial but, in certain situations, may give better outcomes for variations in corneal power^[18-19], surgeon^[12,20], and axial length (e.g., very short or long eyes)^[17,21-24]. Varying lens formula constants according to lens size, to our knowledge, has not been described in the literature.

One limitation to the conclusions drawn from this study is that the statistical differences found in refractive outcomes could be due to differences in axial length rather than the size of the IOL. As there is a correlation between required IOL power and axial length, IOL size may simply be a surrogate marker for axial length. Further studies could look into different sizes of IOLs in eyes with similar axial lengths, for example, patients with bilateral cataracts who have different sized IOLs in each eye. A second limitation is the small number of subjects in the

Table 5 Errors in predicted spherical equivalents (actual postop. spherical equivalent minus predicted spherical equivalent)

Year	Number of cataract operations	Mean error (D)	Standard deviation of error	Errors<0.5 D (%)	Errors<1.0 D (%)	Errors<1.5 D (%)
2008	1409	0.124	0.66	61.4	88.8	96.7
2009	981	0.107	0.613	62.6	90.3	97.4
2010	1265	0.105	0.622	62.2	91.4	97.6
2011	1492	-0.04	0.617	67.3	93.4	97.5
2012	906	0.08	0.584	65.8	94.6	97.4

large IOL group, which limited the inclusion of this group in this study.

ACKNOWLEDGEMENTS

Conflicts of Interest: El-Khayat AR, None; Tesha P, None.

REFERENCES

- Klijn S, Sicam VA, Reus NJ. Long-term changes in intraocular lens position and corneal curvature after cataract surgery and their effect on refraction. *J Cataract Refract Surg* 2016;42(1):35-43.
- Shajari M, Sonntag R, Niermann T, Holland D, Kohnen T, Priglinger S, Mayer WJ. Determining and comparing the effective lens position and refractive outcome of a novel rhexis-fixated lens to established lens designs. *Am J Ophthalmol* 2020;213:62-68.
- Guo HQ, Goncharov A, Dainty C. Intraocular lens implantation position sensitivity as a function of refractive error. *Ophthalmic Physiol Opt* 2012;32(2):117-124.
- Lawu T, Mukai K, Matsushima H, Senoo T. Effects of decentration and tilt on the optical performance of 6 aspheric intraocular lens designs in a model eye. *J Cataract Refract Surg* 2019;45(5):662-668.
- Assia EI, Apple DJ. Side-view analysis of the lens. II. Positioning of intraocular lenses. *Arch Ophthalmol* 1992;110(1):94-97.
- Akreos IOLs : Bausch + Lomb. Accessed on: 2020 Jan 15. <https://www.bausch.com/ecp/our-products/cataract-surgery/lens-systems/akreos-ao>
- Schröder S, Langenbacher A. Relationship between effective lens position and axial position of a thick intraocular lens. *PLoS One* 2018;13(6):e0198824.
- Martinez-Enriquez E, Pérez-Merino P, Durán-Poveda S, Jiménez-Alfaro I, Marcos S. Estimation of intraocular lens position from full crystalline lens geometry: towards a new generation of intraocular lens power calculation formulas. *Sci Rep* 2018;8(1):9829.
- Norrby S. Sources of error in intraocular lens power calculation. *J Cataract Refract Surg* 2008;34(3):368-376.
- Editor CK Senior. Another Step Closer: Lens Constant Customization. Accessed on: 2020 Jan 15. <https://www.reviewofophthalmology.com/article/another-step-closer-lens-constant-customization>
- Ladas JG, Stark WJ. Improving cataract surgery refractive outcomes. *Ophthalmology* 2011;118(9):1699-1700.
- Aristodemou P, Knox Cartwright NE, Sparrow JM, Johnston RL. Intraocular lens formula constant optimization and partial coherence interferometry biometry: refractive outcomes in 8108 eyes after cataract surgery. *J Cataract Refract Surg* 2011;37(1):50-62.
- Ophthalmologists TRC of. Clinical Guidelines. The Royal College of Ophthalmologists. Accessed on: 2020 Jan 19. <https://www.rcophth.ac.uk/standards-publications-research/clinical-guidelines/>
- Overview | Cataracts in adults: management | Guidance | NICE Accessed on: 2020 Jan 15. <https://www.nice.org.uk/guidance/ng77>
- Sheard R. Optimising biometry for best outcomes in cataract surgery. *Eye (Lond)* 2014;28(2):118-125.
- Optimized IOL constants for the ZEISS IOLMaster. Accessed on: 2018 Aug 29. <http://ocusoft.de/ulib/c1.htm>
- Wang L, Koch DD, Hill W, Abulafia A. Pursuing perfection in intraocular lens calculations: III. Criteria for analyzing outcomes. *J Cataract Refract Surg* 2017;43(8):999-1002.
- Eom Y, Kang SY, Song JS, Kim HM. Use of corneal power-specific constants to improve the accuracy of the SRK/T formula. *Ophthalmology* 2013;120(3):477-481.
- Reitblat O, Levy A, Kleinmann G, Lerman TT, Assia EI. Intraocular lens power calculation for eyes with high and low average keratometry readings: comparison between various formulas. *J Cataract Refract Surg* 2017;43(9):1149-1156.
- Langenbacher A, Eppig T, Viestenz A, Seitz B, Müllner G, Schönherr U. Individualization of IOL constants for two hydrophobic intraocular lenses. SRK II, SRK/T, Hoffer-Q, Holladay 1 and Haigis formula. *Ophthalmologie* 2012;109(5):468-473.
- Sanders DR, Retzlaff J, Kraff MC. Comparison of the SRK II formula and other second generation formulas. *J Cataract Refract Surg* 1988;14(2):136-141.
- Zhang JQ, Zou XY, Zheng DY, Chen WR, Sun A, Luo LX. Effect of lens constants optimization on the accuracy of intraocular lens power calculation formulas for highly myopic eyes. *Int J Ophthalmol* 2019;12(6):943-948.
- Gökce SE, Zeiter JH, Weikert MP, Koch DD, Hill W, Wang L. Intraocular lens power calculations in short eyes using 7 formulas. *J Cataract Refract Surg* 2017;43(7):892-897.
- Zheng L, Merriam JC. A method to improve the accuracy of optimized A-constant for IOL calculation formulas. *J Clin Exp Ophthalmol* 2018;9(6):1-3.